

Aptima[®] Multitest Swab Specimen Collection Kit Vaginal Collection Procedure Guide


Swab specimen collection
guide for:

Chlamydia trachomatis (CT)
Neisseria gonorrhoeae (GC)
Trichomonas vaginalis (TV)


Collection for VAGINAL swab specimens: CERVICAL swabs will be rejected

- Not FDA approved for cervical collection
- Incorrect size for cervical collection
- Need mucous removed by cleansing swab prior to cervical collection (not included)


Partially open the swab package and remove the swab. Do not touch the soft tip or lay the swab down. If the soft tip is touched, laid down, or dropped, discard and get a new Aptima Multitest Swab Specimen Collection Kit. Hold the swab, placing thumb and forefinger in the middle of the shaft covering the black score line. Do not hold the shaft below the score line.



Carefully insert the swab into the vagina (NOT cervix) about 2 inches (5 cm) past the introitus and gently rotate the swab for 10 to 30 seconds. Make sure the swab touches the vaginal walls so that moisture is absorbed by the swab. Withdraw the swab without touching the skin.



While holding the swab in hand, unscrew the tube cap. Do not spill the tube contents. **If the tube contents are spilled, discard and replace with a new Aptima Multitest Swab Specimen Collection Kit.**



Immediately place the swab into the transport tube so the black score line is at the top of the tube. Align the score line with the top edge of the tube and carefully break the shaft. The swab will drop to the bottom of the vial. Discard the top portion of the shaft.



Tightly screw the cap onto the tube. When collecting multiple specimens from the same patient, the tube label provides a specimen source field for unique identification for the specimen location.

