



# Marshfield Labs™

A division of Marshfield Clinic

## Targeted Drug Panel, Urine

**\* Quantitations are in units of ng/mL**

|    | Analytes                                    | "Not Detected" | "Present" at "Low Level" | "Present" at "Moderate Level" | "Present" at "High Level" |
|----|---|----------------|--------------------------|-------------------------------|---------------------------|
|    | <b>Amines Drug Class</b>                    |                |                          |                               |                           |
| 1  | Amphetamine                                 | <50            | 50-500                   | 501-2000                      | >2000                     |
| 2  | Methamphetamine                             | <25            | 25-500                   | 501-2000                      | >2000                     |
| 3  | MDA   | <50            | 50-500                   | 501-2000                      | >2000                     |
| 4  | MDMA  | <25            | 25-300                   | 301-600                       | >600                      |
| 5  | MDEA  | <50            | 50-300                   | 301-600                       | >600                      |
| 6  | Methylphenidate                             | <50            | 50-300                   | 301-600                       | >600                      |
| 7  | Ritalinic Acid                              | <50            | 50-300                   | 301-600                       | >600                      |
|    | <b>Benzodiazepine Drug Class</b>            |                |                          |                               |                           |
| 8  | 7-Aminoclonazepam                           | <25            | 25-500                   | 501-3600                      | >3600                     |
| 9  | Lorazepam                                   | <10            | 10-500                   | 501-4363                      | >4363                     |
| 10 | Oxazepam                                    | <10            | 10-500                   | 501-2500                      | >2500                     |
| 11 | Alphahydroxyalprazolam                      | <10            | 10-500                   | 501-2500                      | >2500                     |
| 12 | Nordiazepam                                 | <10            | 10-500                   | 501-1800                      | >1800                     |
| 13 | Temazepam                                   | <10            | 10-500                   | 501-2500                      | >2500                     |
|    | <b>Opioids Drug Class</b>                   |                |                          |                               |                           |
| 14 | Oxymorphone                                 | <10            | 10-500                   | 501-2500                      | >2500                     |
| 15 | Morphine                                    | <25            | 25-1000                  | 1001-4363                     | >4363                     |
| 16 | Hydromorphone                               | <10            | 10-1000                  | 1001-6000                     | >6000                     |
| 17 | Gabapentin                                  | <1000          | 1000-15000               | 15001-50000                   | >50000                    |
| 18 | Codeine                                     | <10            | 10-1000                  | 1001-6000                     | >6000                     |
| 19 | NorTramadol (Tramadol metabolite)           | <15            | 15-500                   | 501-1200                      | >1200                     |
| 20 | 6-Monoacetylmorphine (Heroin metabolite)    | <2             | 2-300                    | 301-800                       | >800                      |
| 21 | Oxycodone                                   | <25            | 25-1000                  | 1001-4363                     | >4363                     |
| 22 | NorOxycodone (Oxycodone metabolite)         | <25            | 25-1000                  | 1001-4363                     | >4363                     |
| 23 | Hydrocodone                                 | <25            | 25-1000                  | 1001-4363                     | >4363                     |
| 24 | Tramadol                                    | <15            | 15-500                   | 501-1200                      | >1200                     |
| 25 | NorBuprenorphine (Buprenorphine metabolite) | <4             | 4-500                    | 501-2500                      | >2500                     |

|    |                                   |     |        |          |       |
|----|-----------------------------------|-----|--------|----------|-------|
| 26 | Buprenorphine                     | <6  | 6-500  | 501-4363 | >4363 |
| 27 | NorFentanyl (Fentanyl metabolite) | <10 | 10-50  | 51-200   | >200  |
| 28 | Fentanyl                          | <4  | 4-500  | 501-1000 | >1000 |
| 30 | EDDP (Methadone Metabolite)       | <50 | 50-500 | 501-1200 | >1200 |
| 31 | Methadone                         | <25 | 25-500 | 501-1200 | >1200 |

JW 9/15/20